

Doyen R100

Bedpan cleaner, mildly alkaline

Application / Dosage

- > Doyen R100 is for commercial applications only.
- > Can be used with the following water hardness: all

Dosage: 1,5 – 3,0 ml / l

- > The dosing is done automatically
- > Use by qualified personnel only

Product Features

- > Cleans intensively
- > Has a high dirt holding capacity
- > Characterized by special economic efficiency

Application

- > Cleaning and disinfection devices, bedpan washers

Important Information

- > Do not mix with other products, do not transfer product into other containers.
- > Tanks, pipes and dosing devices must be rinsed with water before changing products.
- > The operating instructions of the device manufacturer must be observed.

CE marking:  

Doyen R100 complies with European guidelines for medical devices.

If a serious incident occurs with this product, it must be reported to the manufacturer and the relevant national authority.

Doyen R100

Bedpan cleaner, mildly alkaline

Technical Data

- > Colour..... brownish liquid
- > Density..... 1,1 g/ml
- > pH-value 12
- > Storability..... see label imprint
- > Storage conditions.... between 0 and 30°C

Ingredients

Ingredients according to detergent regulations:
5 - 15 % Polycarboxylates
Less than 5 % amphoteric surfactants

- > Chlorine: no
- > Phosphate: no
- > NTA: no

Residual Emptying/ Disposal

- > Empty containers are recyclable and can be disposed of by the end user.
- > Information on the disposal of product residues can be found in the respective safety data sheet.

Sales Units



Can 5 l
Art. no.: 2001485

Hazard information/ Environmental factors



Further information on our products and the EC safety data sheets can be found on :
www.etol.de

