

# Doyen RF90

Cleaner, highly alkaline

## Application / Dosage

- > Doyen RF90 is for commercial applications only.
- > Can be used with the following water hardness: all

**Dosage: 1,0 – 5,0 ml / l**

- > The dosing is done automatically
- > Use by qualified personnel only

## Product Features

- > Cleans intensively
- > Dissolves stubborn soiling
- > Versatile cleaning agent

## Application

- > Particularly suitable for cleaning and disinfection equipment

## Important Information

- > Do not mix with other products, do not transfer product into other containers.
- > Tanks, pipes and dosing devices must be rinsed with water before changing products.
- > The operating instructions of the device manufacturer must be observed.

**Marking: CE MD**

Doyen RF90 complies with European guidelines for medical devices.

If a serious incident occurs with this product, it must be reported to the manufacturer and the relevant national authority.

# Doyen RF90

Cleaner, highly alkaline

## Technical Data

- > Colour..... colourless liquid
- > Density..... 1,3 g/ml
- > pH-value ..... 14
- > Storability..... see label imprint
- > Storage conditions.... between -10 and 30°C

## Ingredients

- Ingredients according to detergent regulations:  
5 - 15 % Phosphonates
- > Chlorine: no
  - > Phosphate: no
  - > NTA: no

## Residual Emptying/ Disposal

- > Empty containers are recyclable and can be disposed of by the end user.
- > Information on the disposal of product residues can be found in the respective safety data sheet.

## Sales Units



Can 5 l  
Art. no.: 2000472



Can 12 kg  
Art. no.: 2001483

## Hazard information/ Environmental factors



Further information on our products and the EC safety data sheets can be found on [www.etol.de](http://www.etol.de)

