

Doyen RS120

Cleaner for cleaning and disinfection appliances, acidic

Application / Dosage

- > **Doyen RS120** is for professional use only
- > Can be used with the following water hardness: all

Dosage: 1.5 - 3.0 ml / l

- > The dosing is done automatically
- > Use by qualified personnel only

Product Features

- > Cleans intensively
- > Reliably dissolves ointment residues
- > Has a high dirt holding capacity
- > Characterized by special economic efficiency

Application

- > Cleaning and disinfection devices, bedpan washers

Important Information

- > Do not mix with other products, do not transfer product into other containers.
- > Tanks, pipes and dosing devices must be rinsed with water before changing products.
- > The operating instructions of the device manufacturer must be observed.

Marking: CE MD

Doyen RS120 complies with European guidelines for medical devices.

If a serious incident occurs with this product, it must be reported to the manufacturer and the relevant national authority.

Doyen RS120

Cleaner for cleaning and disinfection appliances, acidic

Technical Data

- > Colour..... colourless liquid
- > Density..... 1,1 g/ml
- > pH-value 1
- > Storability..... see label imprint
- > Storage conditions.... between 0 and 30°C

Ingredients

Ingredients according to detergent regulations:
15 - 30 % phosphates (phosphoric acid)
Less than 5 % nonionic surfactants, amphoteric surfactants

- > Chlorine: no
- > Phosphate: yes
- > NTA: no

Residual Emptying/ Disposal

- > Empty containers are recyclable and can be disposed of by the end user.
- > Information on the disposal of product residues can be found in the respective safety data sheet.

Sales Units



Can 5 l
Art. no.: 2001486

Hazard information/ Environmental factors



Further information on our products and the EC safety data sheets can be found on www.etol.de

