

etolit Chlortabs

Special cleaner in tablet form, bleaching

Application / Dosage

> **etolit Chlortabs** is suitable for professional only

Dosage: 2 - 5 g / 10 l water (1 tablet = 20 g)

- > The dosing is done manually
- > Add the tablet to the washing tank in addition to the dishwasher detergent
- > For an immersion bath, dissolve in warm water

Product Features

- > Supports machine dishwashing for stubborn coffee and tea stains
- > Suitable for basic cleaning of food colouring deposits
- > Ensures hygienic safety in the rinsing process
- > For machine cleaning or in immersion process

Application

- > With **etolit Chlortabs** the bleaching effect of the dishwasher detergent can be increased in the short term
- > **etolit Chlortabs** can be combined with active-chlorine-free dish cleaners, e.g. when washing the breakfast dishes
- > **etolit Chlortabs** improve the hygienic washing result of the dishwasher
- > Not suitable for silver

Important Information

- > **etolit Chlortabs** is suitable for all food processing companies
- > Do not mix with other products, do not transfer product into other containers
- > Close the lid again after removal!

etolit Chlortabs

Special cleaner in tablet form, bleaching

Technical Data

- > Colour..... White tablets
- > Density..... 1,8 g/ml
- > pH-value 7
- > Storability..... 4 years
- > Storage conditions.... Between -10 and 30°C

Ingredients


Ingredients according to the detergents regulation:
More than 30 % chlorine-based bleaching agent

- > Chlorine: yes
- > Phosphate: no
- > NTA: no

Residual Emptying/ Disposal

- > Empty containers are recyclable and can be disposed of by the end user via DSD.
- > Information on the disposal of product residues can be found in the respective safety data sheet.

Sales Units



Tin 1 x 1 kg
Art.-No.: 2001138
Tin 6 x 1 kg
Art.-No.: 2001132

Hazard information/ Environmental factors



Further information on our products and the EC safety data sheets can be found on :
www.etol.de

