

etolit Freezer Cleaner

Cleaner for freezers and freezer cells, applicable without defrosting

Application / Dosage

- > **etolit Freezer Cleaner** is only for professional use

Dosage: undiluted, up to 10 min exposure time

- > Thick ice crusts should be removed manually before cleaning to ensure the effectiveness of the product

Product Features

- > Ready to use product, immediately applicable
- > No need to clear and defrost the freezer
- > For all alcohol-compatible surfaces
- > Can be used manually or mechanically
- > Highly concentrated and therefore very productive

Application

- > Suitable for manual intermediate cleaning of all alcohol-compatible surfaces such as walls, shelves and doors

Important Information

- > **etolit Freezer Cleaner** is suitable for all food processing companies
- > Do not mix with other products, do not transfer product into other containers
- > Thorough cleaning and disinfection of the deep-freeze rooms can only be achieved when they are defrosted.

etolit Freezer Cleaner

Cleaner for freezers and freezer cells, applicable without defrosting

Technical Data

- > Colour..... Colorless liquid
- > Density..... 1,0 g/ml
- > pH-value 11
- > Storability..... 4 years
- > Storage conditions.... Between -20 and 30°C

Ingredients

Ingredients according to detergent regulations:
Less than 5 % non-ionic surfactants

- > Chlorine: no
- > Phosphate: no
- > NTA: no

Residual Emptying/ Disposal

- > Empty containers are recyclable and can be disposed of by the end user.
- > Information on the disposal of product residues can be found in the respective safety data sheet.

Sales Units



Bottle 4x 1 l
Art.-No.: 2000664



Can 10 l
Art.-No.: 2000520

Hazard information/ Environmental factors



Further information on our products and the EC safety data sheets can be found on :
www.etol.de

