

# etolmat 812

Dosing device for liquid detergents, conductive-controlled device for automatic dosing of etolit cleaners

## Device description

- > Dosing device with conductivity control for automatic input of liquid detergents
- > Self-priming hose pump with spring-loaded rotor rollers and snap-on lid for easy and fast pump hose change
- > Empty display for the cleaner container using conductivity
- > 1 LED (display standby: LED 1 green continuous light, dosing: LED 1 green flashing, deficiency/fault: LED 1 red continuous light)

## Delivery/ Connection material

- 1 Dosing device
- 4 m Connection cable 2 x 0.75 mm<sup>2</sup>
- 4 m Cable for electrode 2 x 0.5 mm<sup>2</sup>
- 1 Sensor 1 pole
- 1 Alkaliean inlet 90°
- 4 m Norprene hose natural ID 6 x AD 12
- 1 Foot valve 12mm
- 1 Hose holder 4.8x6 mm blue
- 1 Hose clamps
- 1 Wall fixing material


## Technical data

- > Voltage..... 230 V, 50 Hz, 16 VA
- > Protection class.....IP 65
- > Dosing area..... 2-75 ml/min
- > Pressure..... Max. 2 bar
- > Suction height.....Max. 1,7m
- > Dimensions..... 65 x 135 x 85 mm

## Application




## Product image



# etolmat 812

## Electrical connection sketch



Datum		21.11.17		etol Oppenau		Projektname		Anlage	
Bearb.				Hygiene		Dosiertechnik		Ort	
Geänd.		04.12.17				etolmat 812		+	
Zustand		Änderung		Datum		Projektbeschreibung		Blatt	
								1	
								von	
								1	

# etolmat 812

## Spare parts

<u>Article</u>	<u>Description</u>
2000524	Console
2000739	Hose holder PP cleaner 6x3 mm blue
2000742	Norprene hose natural ID 6 x AD 12
2002217	Foot valve 12mm
2002230	Electrode single-pole
2002247	Alkaline entry 8mm 90° stainless steel
2003152	Rotor cover 9911 blue
2003173	Pump hose PS 140 No